


# **Cleaning Instructions**

## **Fine Needle Aspirator Handle**



---

**10cc Fine Needle Handle**

Item No. WJ20-0010

**20cc Fine Needle Handle**

Item No. WJ20-0020

---

The F.N.A. handle can be sterilized with ethylene oxide gas or with heated steam at temperature of 121°C for 30 minutes max with 30 minutes drying time.

---



*Manufacturer Update: Limited Autoclavable*

**800-428-1610 • (317) 545-6196**  
**[www.anthonyproducts.com](http://www.anthonyproducts.com)**