

Stainless Steel Vent Tubes MRI Safety

MRI Safety

Our vent tubes with stainless steel or titanium, with or without the stainless wire and other vent tubes (plastic or silicone) with stainless steel wire are MRI Conditional under the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Stainless Steel Vent Tube and or Stainless Steel Wire. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.

Note: Summit Medical vent tubes other than stainless steel or titanium without the stainless steel wire are not affected by MRI's.

Anthony Products, Inc.
7740 Records Street
Indianapolis, IN 46226

Phone 800-428-1610
Fax 317-543-3286 .
Web www.anthonyproducts.com

